



MAMMOGRAPHY

Reading Mammography on Intuitive Imaging's ImageQube Web-based PACS or its Rational Imaging PACS offers both an affordable Mammo solution, as well as a multi-modality workstation.

In a user-friendly environment, Intuitive's Mammography application presents smart templates for both conventional screening and diagnostic exams, and offers the capability of displaying CAD overlays. The structured screening process is effortless, requiring minimal user interaction. The entire Mammo reading process can be handled by the keyboard or the computer mouse.



FEATURES:

- » Intuitive MG screening techniques
- » Flexible Diagnostic view
- » MG Shortcut Toolbar and Thumbnails on each display monitor
- » Anatomic (R, L) and MG view matching (MLO, CC, Film, Stacked or Comparison) hanging protocols
- » Keyboard or mouse controllable smart interface
- » Computer-Aided Detection display capabilities
- » "True Res" magnifying glass and Quadrant View capabilities
- » Robust anatomic annotation features (saved with exam)
- » Zoom, Pan, Rotate, Sort, Magnify and other functionalities
- » FDA-compliant MG printing
- » Auto-Prefetch and Auto-Load comparison exams
- » User-customizable tools and hanging protocols
- » Diagnostic reports, key images and annotations are automatically available to referring physicians through the ImageQube