

PERSONAL IMAGEQUBE

Targeted at radiologists or small imaging centers needing to interpret only low volumes of exams, Intuitive Imaging is pleased to offer the Personal ImageQube.

A single-drive system (no RAID) complete with one Diagnostic license and 5 Clinical licenses, the Personal ImageQube provides PACS technology at an affordable price.



Referring Office

Imaging Center Office

Remote Access

Radiology Department

The Personal ImageQube:

- » Maintains FDA 510(k) clearance
- » Adheres to all imaging standards (HIPAA, DICOM, HL7 and IHE)
- » Receives images and transcribed reports from any standards-based PACS, RIS or imaging acquisition device
- » Runs on any operating system platform (Windows, Macintosh, Solaris or Linux)
- » Provides access to medical images, annotated print pages, transcribed reports or Key Image Summaries
- » Provides data encryption without expensive VPN requirements
- » Offers robust reading functionalities